



Sonic Accessories help hearing instruments do even more

# Add the Wow of Wireless



We envision a world where all people can enjoy the sounds that enrich everyday life. With Sonic, everyday sounds better.



## Optional- or optimal

Wireless technology is all around us, making life more convenient than ever. But did you know you can put it to work via your hearing instruments? Sonic wireless accessories are add-on products that make everyday sounds so accessible, you won't want to listen any other way.



© 2015 Sonic Innovations, Inc. All rights reserved.  
Bluetooth is a registered trademark owned by Bluetooth SIG, Inc, USA.

Apple, the Apple logo, and iPhone are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc. Android and Google Play are trademarks of Google Inc.



# It All Starts with SoundGate 3

## The wireless hub for your hearing instruments

The heart of our wireless system, **SoundGate 3** is the digital “go-between” connecting external devices with your hearing instruments. It’s a small, portable accessory that turns your hearing aids into personal “speakers” for the TV, your phone, or in assistive listening environments.

## Built for Bluetooth®

For example, **SoundGate 3** allows you to connect with any device that’s already Bluetooth® enabled – like your



mobile phone or a digital audio player. Just click a button on SoundGate 3 to transmit the audio through your hearing aids. You keep the SoundGate 3 device in your pocket or wear it around your neck so its microphone can pick up your voice for hands-free phone calls.

## In the loop for telecoil systems

SoundGate 3 is also equipped to take advantage of assistive listening systems. If you visit a theater, classroom, place of worship, or other location that features a hearing loop/telecoil symbol, SoundGate 3 can connect your hearing instruments with the audio there.

## The convenience of remote control

SoundGate 3 offers extra convenience even when you’re not using Bluetooth® or another Sonic accessory because it can function as a remote control for your hearing instruments. Just change programs and volume levels with a discrete click of a button instead of reaching for the controls on your hearing aids themselves.

## Savor every show

With the **TV Adapter 2**, prime time plays right to you. The TV Adapter 2 is an optional accessory that physically connects to your television then wirelessly transmits the audio to SoundGate 3. When you want to watch a show, the audio streams directly to your Sonic hearing devices. It’s like having sports, news, shows and movies play right in your ears.

You control the volume level via SoundGate 3 – turn it up or down to your personal liking. The TV Adapter 2 makes enjoying your favorite programming so convenient, you’ll want an adapter for every TV you own.

## Never miss a call

Answering a call from a mobile phone is already easy with the SoundGate 3 device. But what about your landline phone? That’s where the optional **Phone Adapter 2** comes in. Connect the Phone Adapter 2 to a landline, and it will work with SoundGate 3 to alert you to an incoming call, let you answer remotely, and hear the conversation right through your hearing aids. You can adjust the volume as needed, too. Talking on the phone has never been easier or more clear.

## Connect with more

### SoundGate Mic

With the **SoundGate Mic**, you can close the gap when listening to a distant speaker. It’s a great way to better hear an instructor at the front of the room – or a friend at the far end of the table. Simply clip the SoundGate Mic to the speaker’s clothing or ask them to keep it nearby. Their voice is transmitted via SoundGate 3 right into your hearing instruments.



## SoundLink App for iPhone® and Android™ smartphones

The **SoundLink App**<sup>11</sup> along with SoundGate 3 and all Sonic wireless hearing instruments let you change volume, programs, and listening sources right on your iPhone or Android smartphone. Download the free SoundLink App and let the fun begin. Whether you want to take calls, stream music, or listen to shows, the SoundLink App lets you control how you hear it. Discreet, convenient, and comfortable hearing right at your fingertips.

## RC-N Remote Control

Allows easy operation of program and volume changes. With lower power consumption and easier-to-feel controls, the newly designed **RC-N Remote Control** offers easy, discreet program and volume changes. People with dexterity issues especially love the convenience of this remote – no more fiddling with tiny buttons or dials on hearing instruments.

<sup>11</sup>SoundLink App for iPhone runs on iPhone 6 Plus, iPhone 6, iPhone 5s, iPhone 5c, iPhone 4s, and iPhone 4, and requires Sonic wireless hearing aids and SoundGate 3 with firmware 3.1 or later. SoundLink App for Android smartphones requires Sonic wireless hearing aids and SoundGate 3 with firmware 3.1 or later.